

Multiple-Chamber-Vacuum-Splints

- Division into separate chambers avoids clumping of the granules and assures the immediate initial care
- Compared to the use of pneumatic splints, the immobilisation is pressure-free and therefore avoids pain for the patient.
- Variably positioned Velcro strips (problem-free replacement if lost)
- Comfortable access to heel and elbow
- Well shaped instep
- Full arm splint can also be used as also lower leg splint
- X-ray penetrable, and MRT and CT compatible
- Optional leg immobilisation strap
- SCHNITZLER valve with 7-year manufacturer's guarantee
- Available as a set (lower arm, full arm, leg splint, repair kit, suction pump)

